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Referências

1. Sugrue M et al., World Union of wound Healing Societies (WUWHS) Consensus Document. Closed surgical incision management: understanding the role of NPWT. Wounds International, 2016
2. Jenks et al., (2014) Clinical and economic burden of surgical site infection (SSI) and predicted financial consequences of elimination of SSI from an English hospital, *Journal of Hospital Infection*; **86**, 24-33
3. Badia JM, Casey AL, Petrosillo N, Hudson PM, Mitchel SA. (2017). Impact of surgical site infection on healthcare costs and patient outcomes: a systematic review in six countries. *Journal of Hospital Infection*. 1-15
4. Kirkland, K., Briggs, J., Trivette, S., Wilkinson, W. and Sexton, D. (1999). The impact of surgical site infections in the 1990s: attributable mortality, excess length of hospitalization, and extra costs. *Infection Control and Hospital Epidemiology*, 20(11), pp.725-730
5. Andersson, A., Bergh, I., Karlsson, J. and Nilsson, K. (2010). Patients' experience of acquiring a deep surgical site infection: an interview study. *American Journal of Infection Control*, 38;711-717
6. World Health Organization. (2016) Global guidelines for the prevention of surgical site infection. S.I: World Health Organization, pp.145-148
7. Lumb, H. Data on file report. 1102010. april 2011.
8. Loveluck, J., Copeland, T., Hill, J., Hunt, A., and Martin, R. (2016) Biomechanical modelling of the forces applied to closed incisions during single-use negative pressure wound therapy. *ePlasty*. 16, pp.183-195
9. Malmjsjo, M., Huddleston, E., and Martin, R. (2014) Biological effects of a disposable, canisterless negative pressure wound therapy system. *ePlasty* 10. Canonico S et al., (2012) Therapeutic possibilities of portable NPWT. Initial multidisciplinary observation with the negative pressure device. *Acta Vulnol*. 11. Pp.57-66
11. Pellino et al. (2015) Prophylactic negative pressure wound therapy in colorectal surgery. *Effects of surgical site events: current status and call to action. Updates in Surgery*;67(3):pp235-245
12. Hyldeg, N., Birke-Sorensen, H., Kruse, M., Vinter, C., Joergensen, J., Sorsensen, J., Mogensen, O., Lamont, R. and Bille, C. (2016). Meta-analysis of negative-pressure wound therapy for closed surgical incisions. *British Journal of Surgery*, 103(5), pp.477-486.
13. DS.16.352.R Measurement of NPWT device air leak tolerance – PICO vs PICO September II 2016
14. DS.16.351.R Comparison of noise levels generated by the PICO and PICO II devices in operation. September 2016
15. WMP.11446. UEF/R3 Project Fairbanks Human Factors Summary Report Issue 5. G. Walker, May 2017
16. Hurd, T., Gregory, H., Jones, A. and Brown, S. (2009) A multicentre in-market evaluation of Allevyn Gentle Border; wounds UK, Vol 5, No 3
17. Data on file reference DS/11/037/R2 - In-vitro wound model testing of PICO at a moderate exudate flow rate; Sarah Roberts, March 2011, (in-vitro)
18. Data on file reference DS/11/057/R2 - In-vitro wound model testing of PICO at a low exudate flow rate; Sarah Roberts, April 2011, (in-vitro)
19. Data on file. DS/17/253/R version2. Project Opal PICO 7 system stability testing , initial time point. November 2017
20. Data on file. DS/17/669/R. Negative pressure transmission through dresign port under a compressive force. November 2017
21. Galiano, R., et al., (2014) The effects of a single use canister-free negative pressure wound therapy (NPWT) system* on the prevention of postsurgical wound complications in patients undergoing bilateral breast reduction surgery. Poster presented at The British Association of Aesthetic Plastic Surgeons (BAAPS).
22. Strugala, V. and Martin, R. (2017) Meta-analysis of comparative trails evaluating a prophylactic single-use negative pressure wound therapy system for the prevention of surgical site complications. *Surgical Infections*; DOI: 10.1089/sur.2017.156
23. Tanner et al., (2001), Surveillance of surgical sites in primary care, *Nursing Times*, Vol 107 No 9.
24. Karlakki et al., (2016) Incisional negative pressure wound therapy dressings (iNPWTd) in routing primary hip and knee arthroplasties. *Bone Joint Res* 2016; 5:328-337
25. Bullough, L., Burns, S., Timmons, J., Truman, P. and Megginson, S. (2015) Reducing c-section wound complications. *The Clinical Services Journal*, April, pp.2-6
26. O' Leary, D., Peirce, C., Anglim, B., Carter, M., Hickey, K. Coffey, J. (2016) Prophylactic negative pressure dressing use in closed laparotomy wounds following abdominal operations. *Annals of Surgery*, 265(6), pp.1082-1086
27. Rodden, D and Taylor, A. NPWT: Incision management in high risk cardiothoracic patients – reducing surgical site infection and length of stay. Wounds UK Conference, Harrogate, November 2015
28. Searle R J and Myers D. A survey of caesarean section surgical site infection with PICO single-use negative pressure wound therapy system in high risk patients in England and Ireland. *Journal of hospital infection*. Vol 97, pages 122-124, 2017.
29. Nherera, L, Truman P. and Karlakki, S. (2017) Cost-effectiveness analysis of single-use negative pressure wound therapy dressings (sNPWT) to routine primary hip and knee replacements. *Wound Repair and Regeneration*, (1) pp.1-30
30. Karlakki, S., Brem, M., Giannini, S., Khanduka, V., Stannard, J. and Martin, R. (2013) Negative pressure wound therapy management of the surgical incision in orthopaedic surgery: A review of evidence and mechanisms for an emerging indication. *Bone and Joint Research*, [online] 2(12), pp.276-284. Available at: <http://bjr.boneandjoint.org.uk/content/2/276.long> [Accessed 4 Oct. 2017]
31. Tanner et al., (2009) Post-discharge surveillance to identify colorectal surgical site infection rates and related costs. *Journal of Hospital Infection*. (72), pp.243-250
32. Data on file report DS/14/066/R. Comparison of pressure transmission through RENASYS drape and PICO Soft port dressings. May 2014.

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©Abril 2019 Smith & Nephew